



RNAse Xterminator Spray

#GB43.500S

(For Research only)

Product:	Solution for the removal of RNAse, Enzymes, and DNA from laboratory surfaces
Quantity:	500ml
Appearance:	aqueous, transparent
Storage:	room temperature

Usage:

RNAse Xterminator Spray is a secondary alcohol based solution containing anionic surfactants for the inactivation and removal of RNAses (and other enzymes) as well as DNA from laboratory surfaces. The solution can be used to decontaminate bench tops, labware, pipettors, and other equipment. It is effective on metal, glass and plastic surfaces. Available as a ready-to-use solution, simply spray on the contaminated area and wipe away from the surface using ultrapure water.