

# Mobicol "F" with fixed outlet plug and 2 different screw caps, without filters




## Product Information Sheet # M1050



### SUMMARY

shipped at RT; store at RT

For research use only



Syringe connected to Luer-lock cap


Mobicol "F" is a 1 ml column with particular handling advantages.

Mobicol "F"

- is compatible to the laboratory standard,
- can be centrifuged in minifuges,
- has a Luer-lock adapter,
- smallest and largest volumes can be treated easily.



Filters of the pore sizes 10  $\mu\text{m}$ , 35  $\mu\text{m}$  or 90  $\mu\text{m}$  are available for Mobicols. The filters can be inserted easily with the help of a plastic tool which is delivered with the columns. The upper large filter fits tightly over the whole length of the column body.

Mobicols can be closed tightly with a screw-on cap and the reversed side of the snap-off plug at the outlet. The volume is 700  $\mu\text{l}$ . Mobicols can be autoclaved at 120  $^{\circ}\text{C}$  (with inserted filter at 110  $^{\circ}\text{C}$  for 10 minutes).



Adapter M3001: tubing to Mobicol

### Mobicol "F"




Mobicols can be centrifuged (fit into 1.5 ml or 2.2 ml tubes)

# Mobicol "F" with fixed outlet plug and 2 different screw caps, without filters



## Product Information Sheet # M1050



The Mobicol "F" M1050 is supplied with two screw caps (Luer-lock cap and screw cap) and a snap-off plug for the outlet. Filters with different pore sizes (10  $\mu\text{m}$ , 35  $\mu\text{m}$  and 90  $\mu\text{m}$ ) are separately available, so you can use the filter of your choice. With this features the Mobicol "F" can be used for a variety of applications.

### Order Information, Shipping and Storage

| Order#   | Product   | Quantity |
|----------|---|----------|
| M1050    | Mobicol "F" with fixed outlet plug and 2 different screw caps, without filters                  | 50       |
| M1053    | Mobicol "F" with fixed outlet plug and screw cap, without filters                               | 50       |
| M105035F | Mobicol "F" with fixed outlet plug, inserted 35 $\mu\text{m}$ filter and 2 different screw caps | 50       |
| M2210    | Filter (large) 10 $\mu\text{m}$ pore size, for Mobicol M1002, M1003, M1050 & M1053              | 50       |
| M523515  | Filter (large) 35 $\mu\text{m}$ pore size, for Mobicol M1002, M1003, M1050 & M1053              | 50       |
| M2290    | Filter (large) 90 $\mu\text{m}$ pore size, for Mobicol M1002, M1003, M1050 & M1053              | 50       |

shipped at RT; store at RT

### Column Accessories

| Order# | Product  | Quantity |
|--------|--|----------|
| M3001  | Tubing-to-Luer-lock adapter (material outside: metal; inside: glass and PTFE; with 1 m PTFE tubing) for Mobicols or laboratory columns | 1        |
| M3002  | Luer adapter, (M+F) 10 each  | 20       |
| M3003  | Luer adapter for the top-cap (M)   | 20       |
| M3004  | Luer adapter for outlet (F)  | 20       |
| M3006  | Stopcocks, one-way, Luer-lock  | 4        |
| M3007  | Luer-lock adapter (M/M)  | 1        |

shipped at RT; store at RT